




by **NEWTOP**



**DIKTE
VERGELIJK
EN TOOLS**

STANDAARD RX- EN VOORRAADGLAZEN

65 Ø	+ 2.00 dpt	+ 4.00 dpt	+ 6.00 dpt	+ 8.00 dpt
1.50 UNV	0,8 mm CT= 3,0 mm	0,8 mm CT= 5,1 mm	0,8 mm CT= 7,28 mm	0,8 mm CT= 9,51 mm
1.60 UNV	0,8 mm CT= 2,6 mm	0,8 mm CT= 4,4 mm	0,8 mm CT= 6,15 mm	0,8 mm CT= 7,98 mm
1.67 UNV	0,8 mm CT= 2,45 mm	0,8 mm CT= 3,98 mm	0,8 mm CT= 5,58 mm	0,8 mm CT= 7,21 mm
1.74 UNV	1,2 mm CT= 2,62 mm	0,9 mm CT= 3,75 mm	0,9 mm CT= 5,24 mm	0,9 mm CT= 6,71 mm














INDIVIDUELE RX-GLAZEN

	65 Ø	+ 2.00 dpt	+ 4.00 dpt	+ 6.00 dpt	+ 8.00 dpt
1.50 PCS		0,8 mm CT= 2,89 mm	0,8 mm CT= 4,89 mm	0,7 mm CT= 6,74 mm	0,7 mm CT= 8,83 mm
1.60 PCS		0,8 mm CT= 2,57 mm	0,8 mm CT= 4,04 mm	0,7 mm CT= 5,67 mm	0,7 mm CT= 7,46 mm
1.67 PCS		0,9 mm CT= 2,37 mm	0,8 mm CT= 3,62 mm	0,8 mm CT= 5,19 mm	0,8 mm CT= 6,83 mm
1.74 PCS		1,4 mm CT= 2,77 mm	0,8 mm CT= 3,38 mm	0,8 mm CT= 4,83 mm	0,8 mm CT= 6,43 mm


STANDAARD RX- EN VOORRAADGLAZEN

70 Ø	- 2.00 dpt	- 4.00 dpt	- 6.00 dpt	- 8.00 dpt
1.50	<p>4,67 mm CT= 2,0 mm</p>	<p>7,15 mm CT= 1,90 mm</p>	<p>10,06 mm CT= 1,80 mm</p>	<p>12,70 mm CT= 1,80 mm</p>
1.60	<p>3,7 mm CT= 1,6 mm</p>	<p>5,92 mm CT= 1,6 mm</p>	<p>8,12 mm CT= 1,6 mm</p>	<p>10,43 mm CT= 1,6 mm</p>
1.67	<p>3,44 mm CT= 1,5 mm</p>	<p>5,35 mm CT= 1,5 mm</p>	<p>7,34 mm CT= 1,5 mm</p>	<p>9,32 mm CT= 1,5 mm</p>
1.74	<p>3,3 mm CT= 1,6 mm</p>	<p>5,1 mm CT= 1,6 mm</p>	<p>6,99 mm CT= 1,6 mm</p>	<p>8,69 mm CT= 1,6 mm</p>

INDIVIDUELE RX-GLAZEN


70 Ø	- 2.00 dpt	- 4.00 dpt	- 6.00 dpt	- 8.00 dpt
1.50 UNV	 <p>4,55 mm</p> <p>CT= 2,0 mm</p>	 <p>6,87 mm</p> <p>CT= 1,90 mm</p>	 <p>9,95 mm</p> <p>CT= 1,90 mm</p>	 <p>11,92 mm</p> <p>CT= 1,80 mm</p>
1.60 UNV	 <p>3,40 mm</p> <p>CT= 1,4 mm</p>	 <p>5,40 mm</p> <p>CT= 1,4 mm</p>	 <p>7,35 mm</p> <p>CT= 1,4 mm</p>	 <p>9,45 mm</p> <p>CT= 1,35 mm</p>
1.67 UNV	 <p>3,16 mm</p> <p>CT= 1,35 mm</p>	 <p>4,87 mm</p> <p>CT= 1,35 mm</p>	 <p>6,73 mm</p> <p>CT= 1,35 mm</p>	 <p>8,44 mm</p> <p>CT= 1,35 mm</p>
1.74 UNV	 <p>3,04 mm</p> <p>CT= 1,35 mm</p>	 <p>4,61 mm</p> <p>CT= 1,35 mm</p>	 <p>6,55 mm</p> <p>CT= 1,35 mm</p>	 <p>8,24 mm</p> <p>CT= 1,35 mm</p>


Design	Gravering	Design symbolen		Extra symbool		Additie	Index
Compass	o	Cm	14 / 15 / 16 / 17 / 18	-	-	-	1,50/1,53/1,59/1,60/1,67/1,74
Effecto	o	E	14 / 15 / 16 / 17 / 18	-	p	ADD	1,50/1,53/1,59/1,60/1,67/1,74
Velveto	o	VL	14 / 15 / 16 / 17 / 18	-	p	ADD	1,50/1,53/1,56/1,59/1,60/1,67/1,74
Full Screen	o	MF	14 / 15 / 16 / 17 / 18	-	p	ADD	1,50/1,53/1,56/1,59/1,60/1,67/1,74
MultiFit H	o	MH	14 / 15 / 16 / 17 / 18	-	p	ADD	1,50/1,53/1,56/1,59/1,60/1,67/1,74
MultiFit S	o	MS	14 / 15 / 16 / 17 / 18	-	p	ADD	1,50/1,53/1,56/1,59/1,60/1,67/1,74
Natura	o	NA	14 / 16 / 18 / 20	-	-	ADD	1,50/1,53/1,56/1,59/1,60/1,67/1,74
Camber Office	o	c	13 / 2 / 4 / 6	4/8	-	ADD	1,50/1,53/1,59/1,60/1,67/1,74
Office	o	13 / 2 / 4 / 6	M	4/8	p	ADD	1,50/1,53/1,56/1,59/1,60/1,67/1,74
Office VIDE	o	V	4 / 8	-	-	ADD	1,50/1,60
EyeChi	o	E	-	-	p	25 / 50 / 75 / 1,0	1,50/1,53/1,56/1,59/1,60/1,67/1,74
No Tense III	o	nT	-	-	p	25 / 50 / 75 / 1,0	1,50/1,53/1,56/1,59/1,60/1,67/1,76
Okos	o	oK	-	-	p	-	1,50/1,53/1,56/1,59/1,60/1,67/1,75
Bifo Digital	o	Bi	-	-	-	ADD	1,50/1,60/1,67
Bifo Round 30	o	B	-	-	-	ADD	1,50/1,53/1,59/1,60/1,67
PCS	o	PCS	-	-	-	-	1,50/1,53/1,56/1,59/1,60/1,67/1,74



Hoogte

Uitlijning





**HVA omreken tabel
op 14 mm**

MINUS	0 mm	PLUS
-3.13	3.00	+2.88
-3.68	3.50	+3.34
-4.24	4.00	+3.79
-4.80	4.50	+4.23
-5.38	5.00	+4.67
-5.67	5.25	+4.89
-5.96	5.50	+5.11
-6.25	5.75	+5.32
-6.55	6.00	+5.54
-6.85	6.25	+5.75
-7.15	6.50	+5.96
-7.45	6.75	+6.17
-7.76	7.00	+6.38
-8.07	7.25	+6.58
-8.38	7.50	+6.79
-8.70	7.75	+6.99
-9.01	8.00	+7.19
-9.65	8.50	+7.60
-10.30	9.00	+7.99
-10.96	9.50	+8.39
-11.63	10.00	+8.77
-12.31	10.50	+9.15
-13.00	11.00	+9.53

